



Product Datasheet

ANXA1 Antibody [SAIC-13B-19], Rabbit IgG (orb1672072)


| | |
|------------------------|--|
| Description | ANXA1 Antibody [SAIC-13B-19], Rabbit IgG |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Immunogen | Peptide "AAYLQETGKPLDETLK" derived from Annexin A1 conjugated to KLH. |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | batch dependent |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Isotype | IgG kappa |
| Clonality | Recombinant |
| Uniprot ID | P04083 |
| Expiration Date | 12 months from date of receipt. |



Western Blot using anti-ANXA1 antibody S...



Flow cytometry using the Anti-ANXA1 anti...



Immunofluorescence staining of HeLa cell...