



Product Datasheet

CALM1P2 Antibody (orb1672076)

Description	CALM1P2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide "TFIPVTSAK" derived from CD44, a hyaluronic acid receptor mediating cell-cell and cell-matrix interactions , conjugated to KLH.
Target	CALM1P2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P16070
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Flow cytometry using the Anti-CD44 (Indi...