



Product Datasheet

SNAP-25 Antibody [SP12], Mouse IgG1 (orb1672080)

Description	SNAP-25 Antibody [SP12], Mouse IgG1
Reactivity	Bovine, Feline, Guinea pig, Hamster, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rin usin...



Western Blot using anti-SNAP25 antibody

...