

Product Datasheet

CD33 Biosimilar Antibody (orb1672100)

Description	CD33 Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Gemtuzumab was generated by humanization of the murine P67.6 clone, raised against CD33-transfected FMY9S5 cells of human origin. CD33 is a transmembrane receptor expressed on predominantly on cells of the myeloid lineage, but it can also be found on some lymphoid cells.
Target	CD33
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG4 kappa
Clonality	Recombinant
Uniprot ID	P20138
Expiration Date	12 months from date of receipt.