

## Product Datasheet

# carcinoembryonic antigen (CEA) Antibody [Arcitumomab], Rabbit IgG - Research Grade Biosimilar (orb1672103)

<b>Description</b>	carcinoembryonic antigen (CEA) Antibody [Arcitumomab], Rabbit IgG - Research Grade Biosimilar
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IHC-P
<b>Immunogen</b>	Carcinoembryonic antigen (CEA) describes a set of highly related glycoproteins involved in cell adhesion normally produced in gastrointestinal tissue during fetal development, but only at very low levels in adults.
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q13984</a>
<b>Expiration Date</b>	12 months from date of receipt.