



## Product Datasheet

### C3 Biosimilar Antibody (orb1672108)

<b>Description</b>	C3 Biosimilar Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, FC
<b>Immunogen</b>	The epidermal growth factor receptor is the cell-surface receptor for members of the epidermal growth factor family (EGF-family) of extracellular protein ligands.
<b>Target</b>	C3
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1 kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P00533</b>
<b>Expiration Date</b>	12 months from date of receipt.



Binding curve of anti-EGFR antibody to E...



Immunohistochemical staining of rat stom...