



Product Datasheet

CD7 Antibody (orb1672110)

Description	CD7 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC
Target	CD7
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P09564
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed Jurkat cells.



Western Blot using anti-CD7 antibody 3A1.