



Product Datasheet

CR1 Antibody (orb1672121)

Description	CR1 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, COSTIM, FC, IHC-Fr, IP, WB
Immunogen	Purified Mouse CR1.
Target	CR1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P19070, P17927
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat stom...



Flow-cytometry using anti-complement rec...