


## Product Datasheet

# Flavivirus group antigen Antibody [D1-4G2-4-15 (4G2)], Rabbit IgG (orb1672124)

<b>Description</b>	Flavivirus group antigen Antibody [D1-4G2-4-15 (4G2)], Rabbit IgG
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, IF, IHC-P, WB
<b>Immunogen</b>	Dengue Virus type 2 antigens.
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.



Detection of Zika virus by immunofluores...