


Product Datasheet

Flavivirus group antigen Antibody [D1-4G2-4-15 (4G2)], Mouse IgG2a (orb1672125)

Description	Flavivirus group antigen Antibody [D1-4G2-4-15 (4G2)], Mouse IgG2a
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-P, WB
Immunogen	Dengue Virus type 2 antigens.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



ELISA on Denge Virus-Like Particles usin...