


Product Datasheet

AFP Antibody (orb1672138)


Description	AFP Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	AF5 was prepared by immunizing female BALB/c mice with human amniotic fluid containing AFP, from a pregnancy affected by a neural tube defect, diluted and emulsified in Freund's complete adjuvant.
Target	AFP
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P02771
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-Alpha-feto...



Western Blot using anti-Alpha-fetoprotei...



Immunofluorescence staining of fixed HepG...