



## Product Datasheet

### FAS Biosimilar Antibody (orb1672140)


<b>Description</b>	FAS Biosimilar Antibody
<b>Reactivity</b>	Human, Primate
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, FC, IF, IHC, WB
<b>Immunogen</b>	Basiliximab was prepared by immunizing BALB/c mice with CTC 16 cells.
<b>Target</b>	FAS
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P01589</b>
<b>Expiration Date</b>	12 months from date of receipt.



Flow-cytometry using the anti-CD25 (IL2R...



Flow-cytometry using anti-CD25 antibody ...



Flow-cytometry using anti-CD25 antibody ...