



Product Datasheet

FAS Biosimilar Antibody (orb1672141)


Description	FAS Biosimilar Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IF, IHC, WB
Immunogen	Basiliximab was prepared by immunizing BALB/c mice with CTC 16 cells.
Target	FAS
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P01589
Expiration Date	12 months from date of receipt.



Western Blot using anti-IL2R alpha (CD25...



Immunohistochemical staining of rat sple...



Immunofluorescent staining of rat thymus...