



Product Datasheet

FAS Biosimilar Antibody (orb1672142)


Description	FAS Biosimilar Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	Cultured T-cells derived from a patient with Mycosis fungoides.
Target	FAS
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P01589
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-CD25 (IL2R...



Flow-cytometry using anti-IL-2R antibody...



Flow-cytometry using anti-IL-2R antibody...