
Product Datasheet

**HLA class II Antibody [F3.3], Rabbit IgG
(orb1672144)**

Description

HLA class II Antibody [F3.3], Rabbit IgG

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, WB

Preservatives

PBS with 0.02% Proclin 300.

Concentration

batch dependent

Storage

Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.

Note

For research use only

Isotype


IgG kappa

Clonality


Recombinant

Expiration Date


12 months from date of receipt.



Western Blot using anti-HLA class II his...



Flow-cytometry using the anti-HLA II DR ...



Flow-cytometry using anti-HLA class II a...