

Product Datasheet

EEC1 Antibody (orb1672148)

Description	EEC1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Depletion, IF, IHC, WB
Immunogen	Molt 4, a human T cell line.
Target	EEC1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b kappa
Clonality	Recombinant
Uniprot ID	P61769
Expiration Date	12 months from date of receipt.