

Product Datasheet

Unknown Specificity (Isotype control) Antibody [MOPC-21], Rat IgG1 (orb1672155)

| | |
|----------------------------|---|
| Description | Unknown Specificity (Isotype control) Antibody [MOPC-21], Rat IgG1 |
| Conjugation | Unconjugated |
| Tested Applications | Control |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | batch dependent |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Isotype | IgG1 kappa |
| Clonality | Recombinant |
| Expiration Date | 12 months from date of receipt. |