

Product Datasheet

Unknown Specificity (Isotype control) Antibody [MOPC-21], Mouse IgG3 (orb1672156)

Description	Unknown Specificity (Isotype control) Antibody [MOPC-21], Mouse IgG3
Conjugation	Unconjugated
Tested Applications	Control
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG3 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.