



Product Datasheet

NOTCH1 Antibody (orb1672173)

Description	NOTCH1 Antibody
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	This antibody was selected against a fusion protein consisting of EGF domains 1-12 of murine Notch1 fused to a human Fc domain and binds to to the NRR of mouse Notch1.
Target	NOTCH1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	P46531 , Q01705
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-Notch 1 an...



Immunofluorescence staining of fixed HeLa...