



Product Datasheet

CD3E Antibody (orb1672178)


Description	CD3E Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human CD3.
Target	CD3E
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	P07766
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD3 antibody Y...



Flow-cytometry using the anti-CD3E antib...



Western Blot using anti-CD3E antibody YT...