



Product Datasheet

CD4 Antibody (orb1672183)

Description	CD4 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-Fr
Immunogen	Rat T-cell line (NB2) transfected with the human CD4.
Target	CD4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.



Western Blot using anti-CD4 antibody. Mo...



Immunohistochemical staining of reen usi...