



Product Datasheet

Cd8a Antibody (orb1672186)

| | |
|----------------------------|--|
| Description | Cd8a Antibody |
| Reactivity | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | Depletion, FC, IHC, IHC-Fr, In vivo |
| Immunogen | Murine CD8 (Iyt-2). |
| Target | Cd8a |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | batch dependent |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Isotype | IgG kappa |
| Clonality | Recombinant |
| Uniprot ID | P01731 |
| Expiration Date | 12 months from date of receipt. |



Flow-cytometry using anti-CD8alpha/Lyt-2...



Flow-cytometry using the anti-CD8 alpha/...