


Product Datasheet

CD2 Antibody (orb1672191)

Description	CD2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Depletion, FC
Target	CD2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P06729
Expiration Date	12 months from date of receipt.



Western Blot using anti-CD2 antibody YTH...



Immunofluorescence staining of fixed Mol...