

## Product Datasheet

### IGKV1OR1-1 Antibody (orb1672203)

<b>Description</b>	IGKV1OR1-1 Antibody
<b>Reactivity</b>	Human, Monkey, Mouse, Rabbit, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-Fr, IHC-P, IP, In vivo, WB
<b>Immunogen</b>	Synthetic peptide corresponding to aa 371-380 of human p53.
<b>Target</b>	IGKV1OR1-1
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">P04637</a>
<b>Expiration Date</b>	12 months from date of receipt.