

Product Datasheet

V5 epitope tag Antibody [SV5-P-K], Human Fab fragment (orb1672223)

Description	V5 epitope tag Antibody [SV5-P-K], Human Fab fragment
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IP, RIA, WB
Immunogen	Pk (V5) Epitope Tag (GKPIP NPLLGLDST).
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Tag	His-Tagged
Note	For research use only
Isotype	Fab fragment kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.