

## Product Datasheet

# V5 epitope tag Antibody [SV5-P-K], Mouse Fab fragment (orb1672225)

<b>Description</b>	V5 epitope tag Antibody [SV5-P-K], Mouse Fab fragment
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-Fr, IP, RIA, WB
<b>Immunogen</b>	Pk (V5) Epitope Tag (GKPIP NPLLGLDST).
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Tag</b>	His-Tagged
<b>Note</b>	For research use only
<b>Isotype</b>	Fab fragment kappa
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.