


Product Datasheet

V5 epitope tag Antibody [SV5-P-K], Mouse IgG1 (orb1672226)

Description	V5 epitope tag Antibody [SV5-P-K], Mouse IgG1
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IP, RIA, WB
Immunogen	Pk (V5) Epitope Tag (GKPIPPLLGLDST).
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



ELISA of anti-epitope tag antibodies on ...