



Product Datasheet

CLEC4M Antibody (orb1672248)

Description	CLEC4M Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	CD299
Target	CLEC4M
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9H2X3
Expiration Date	12 months from date of receipt.



Western Blot using anti-DC-SIGNR (CLEC4M...



Immunofluorescence staining of fixed HeLa...