



Product Datasheet

CD34 Antibody (orb1672254)


Description	CD34 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	KG-1a Human myeloblastic leukemic cell line (immature human myeloid cell line) from a patient with nonlymphocytic leukemia.
Target	CD34
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P28906
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD34 antibody ...



Immunohistochemical staining of human pl...



Western Blot using anti-CD34 antibody My...