



Product Datasheet

MS4A1 Biosimilar Antibody (orb1672256)


Description	MS4A1 Biosimilar Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human lymphoblastoid cell line SB.
Target	MS4A1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P11836
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-CD20 resea...



Flow-cytometry using the anti-CD20 resea...



Flow-cytometry using the anti-CD20 resea...