



Product Datasheet

Fcgr2b Antibody (orb1672260)

Description	Fcgr2b Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IF, IHC, IP
Immunogen	Mouse BALB/c Macrophage J774 Cell Line.
Target	Fcgr2b
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b kappa
Clonality	Recombinant
Uniprot ID	P08101, P08508
Expiration Date	12 months from date of receipt.



Blocking of Fc-receptors by anti-CD16/32...



Anti-CD16/32 staining of murine macroph...