

## Product Datasheet

# His tag (C-term) Antibody [3D5], Rat IgG2a (orb1672315)

<b>Description</b>	His tag (C-term) Antibody [3D5], Rat IgG2a
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Purified his tagged Fab fragment.
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a kappa
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.