

Product Datasheet

His tag (C-term) Antibody [3D5], Hamster (Armenian) IgG (orb1672321)

Description	His tag (C-term) Antibody [3D5], Hamster (Armenian) IgG
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified his tagged Fab fragment.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.