

## Product Datasheet

### His tag (C-term) Antibody [3D5], Human IgG3 (orb1672324)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | His tag (C-term) Antibody [3D5],<br>Human IgG3                                   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, WB  |
| <b>Immunogen</b>           | Purified his tagged Fab fragment.  |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.  |
| <b>Concentration</b>       | batch dependent  |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG3 kappa   |
| <b>Clonality</b>           | Recombinant  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |