



## Product Datasheet

### MYC Antibody (orb1672344)

<b>Description</b>	MYC Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC-Fr, IHC-P, IP, WB
<b>Immunogen</b>	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
<b>Target</b>	MYC
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1 kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P01106</b>
<b>Expiration Date</b>	12 months from date of receipt.



Left Panel:  
orb1672344  
(2 ug/ml)  
stainin...



ELISA of anti-  
epitope tag  
antibodies on

...