



Product Datasheet

ZC3H12A Antibody (orb1672348)


Description	ZC3H12A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Anti-ZC3H12A antibody (orb1672348) was raised against a peptide corresponding to 14 amino acids near the carboxyl terminus of human ZC3H12A. The immunogen is located within the last 50 amino acids of ZC3H12A.
Target	ZC3H12A
Preservatives	ZC3H12A Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	ZC3H12A antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 66 kD Observed: 66 kD
Uniprot ID	Q5D1E8
NCBI	Q5D1E7.2
Expiration Date	12 months from date of receipt.



WB
Validation in
Human Lung.
Loading:
15...



Western Blot
Validation in
Mouse
Tissues...



Western Blot
Validation in
Rat Tissues.
...