

Product Datasheet

Recombinant Human Maspin (SerpinB5) (orb1673606)

Description	Recombinant Human Maspin (SerpinB5)
Storage	Store product at -70°C . For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favourable performance, avoid repeated handling and multiple freeze/thaw cycles.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	42.3 kDa
Expiration Date	6 months from date of receipt.