

Product Datasheet

Recombinant Human IL36RN (IL1F5) (orb1673690)

Description	Recombinant Human IL36RN (IL1F5)
Form/Appearance	Lyophilized from a 0.2 µm filtered solution in Tris, NaCl and KCl
Concentration	1.0 EU/µg of recombinant protein as determined by the LAL method
Storage	This protein can be stored in working aliquots at 2° - 8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
Solubility (25°C)	A quick spin of the vial followed by reconstitution in sterile distilled water to a concentration not less than 0.1 mg/mL is recommended. This solution can then be diluted into other buffers.
Expiration Date	6 months from date of receipt.