

Product Datasheet

Recombinant Human IGFBP7 (orb1673767)

| | |
|--------------------------|---|
| Description | Recombinant Human IGFBP7 |
| Form/Appearance | Lyophilized from a 0.2 µm filtered solution in PBS |
| Concentration | 1.0 EU/µg of recombinant protein as determined by the LAL method. |
| Storage | 1.0 EU/µg of recombinant protein as determined by the LAL method. |
| Note | For research use only |
| Purity | > 95% as determined by SDS-PAGE |
| Source | E. coli |
| MW | 26.7 kDa |
| Solubility (25°C) | Reconstitute at 100 µg/mL in sterile PBS |
| Expiration Date | 6 months from date of receipt. |