

## Product Datasheet

# Recombinant Human ERN1 (IRE1) (N-GST tag) (orb1673917)

<b>Description</b>	Recombinant Human ERN1 (IRE1) (N-GST tag)
<b>Form/Appearance</b>	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
<b>Storage</b>	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	70% - 90%
<b>Source</b>	Insect cells
<b>MW</b>	87 kDa
<b>Expiration Date</b>	6 months from date of receipt.