

Product Datasheet

Recombinant Human BMP-3 beta (GDF10) (orb1674095)

Description	Recombinant Human BMP-3 beta (GDF10)
Form/Appearance	Recombinant BMP-3b was lyophilized from a 0.2 µm filtered 1.0 mg/mL 10 mM sodium citrate (pH 3.5) containing 10% glycerol.
Concentration	1.0 EU/µg of recombinant protein as determined by the LAL method.
Storage	Upon reconstitution, this cytokine can be stored in working aliquots at 2° - 8°C for one month, or at -20°C for six months. Avoid repeated freeze/thaw cycle.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	13 kDa
Solubility (25°C)	A quick spin of the vial followed by reconstitution in distilled water to a concentration not less than 0.1 mg/mL. This solution can then be diluted into other buffers.
Expiration Date	6 months from date of receipt.