

## Product Datasheet

### Recombinant Human Activin A (INHBA) (orb1674169)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Recombinant Human Activin A (INHBA)  |
| <b>Form/Appearance</b>   | Activin A is supplied in 20mM Tris-HCl, 1mM EDTA and 50% Glycerol.                                   |
| <b>Concentration</b>     | 1.0 EU/ $\mu$ g of recombinant protein as determined by the LAL method.                              |
| <b>Storage</b>           | Store at 4°C, if the product will be used within 4 weeks. Store at -20°C for longer periods of time. |
| <b>Note</b>              | For research use only  |
| <b>Purity</b>            | > 95% as determined by SDS-PAGE  |
| <b>Source</b>            | E. coli  |
| <b>MW</b>                | 13 kDa   |
| <b>Solubility (25°C)</b> | Not required   |
| <b>Expiration Date</b>   | 6 months from date of receipt.   |