

Product Datasheet

IgG Fraction of anti-Rabbit IgG [H&L] [Goat] (orb1674971)

Description	IgG Fraction of anti-Rabbit IgG [H&L] [Goat]
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	Conjugation, ID, IP
Immunogen	Rabbit IgG whole molecule
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	10.0 mg/ml (by UV absorbance at 280 nm)
Storage	Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.
Note	For research use only
Application notes	ELISA 1:70,000 - 1:300,000 WESTERN BLOT 1:5,000 - 1:25,000 IMMUNOHISTOCHEMISTRY 1:1,000 - 1:5,000 OTHER APPLICATIONS User Optimized
Expiration Date	12 months from date of receipt.