

Product Datasheet

Agarose Immobilized IgG Fraction of anti-Rat IgG [H&L] [Goat] (orb1674976)

Description	Agarose Immobilized IgG Fraction of anti-Rat IgG [H&L] [Goat]
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	Purification
Immunogen	Rat IgG whole molecule
Preservatives	0.02 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.3
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE. Expiration date is one (1) year from date of opening.
Note	For research use only
Application notes	Not specified. Optimal titers for applications should be determined by the researcher.
Expiration Date	12 months from date of receipt.