

Product Datasheet

Agarose Immobilized IgG Fraction of anti-Rabbit IgG (H&L) [Goat] (orb1674977)

| | |
|----------------------------|--|
| Description | Agarose Immobilized IgG Fraction of anti-Rabbit IgG (H&L) [Goat] |
| Species/Host | Goat |
| Reactivity | Rabbit |
| Conjugation | Unconjugated |
| Tested Applications | Purification |
| Immunogen | Rabbit IgG whole molecule |
| Preservatives | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Storage | Store vial at 4° C prior to opening. DO NOT FREEZE. Expiration date is one (1) year from date of opening. |
| Note | For research use only |
| Application notes | Not specified. Optimal titers for applications should be determined by the researcher. |
| Expiration Date | 12 months from date of receipt. |