

Product Datasheet

Fluorescein Conjugated IgG fraction of anti-Mouse IgM [H&L] [Rabbit] (orb1674982)

Description	Fluorescein Conjugated IgG fraction of anti-Mouse IgM [H&L] [Rabbit]
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FACS, FC
Immunogen	Mouse IgM, whole molecule
Preservatives	0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.6
Concentration	10.0 mg/ml (by UV absorbance at 280 nm)
Storage	Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of restoration.
Note	For research use only
Application notes	Not specified. Optimal titers for applications should be determined by the researcher
Expiration Date	12 months from date of receipt.