
Product Datasheet

Peroxidase Conjugated Affinity Purified anti-Donkey IgG [H&L] [Rabbit] (orb1675014)

Description

Peroxidase Conjugated Affinity

	[Rabbit]
Species/Host	Rabbit
Reactivity	Donkey
Conjugation	HRP
Tested Applications	DOT, ELISA, EM, IHC-P, WB
Immunogen	Donkey IgG whole molecule
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	2.0 mg/ml (by UV absorbance at 280 nm)
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
Note	For research use only
Application notes	This product has been assayed against 1.0 µg of Donkey IgG in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:20,000 to 1:100,000 of the reconstitution concentration is suggested for this product. Optimal titers for other applications should be determined by the researcher.
Expiration Date	12 months from date of receipt.