

Product Datasheet

Anti-HUMAN IgM (mu chain) (GOAT) Antibody Peroxidase Conjugated (orb1675024)

Description	Anti-HUMAN IgM (mu chain) (GOAT) Antibody Peroxidase Conjugated
Species/Host	Goat
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Human IgM(mu) heavy chain
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	1.0 mg/ml (by UV absorbance at 280 nm)
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
Note	For research use only
Application notes	ELISA: 1:180,000 / WB: 1:2,000 - 1:20,000 / IHC: 1:500 - 1:2,500 / OTHER: User Optimized
Expiration Date	12 months from date of receipt.