

---

## Product Datasheet

**Cy3 Conjugated Affinity Purified anti-Mouse IgG  
[H&L] [Goat] (orb1675059)**

**Description**

Cy3 Conjugated Affinity Purified

<b>Species/Host</b>	Goat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Cy3.29
<b>Tested Applications</b>	FACS, FC
<b>Immunogen</b>	Mouse IgG, whole molecule
<b>Preservatives</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Concentration</b>	1.0 mg/ml (by UV absorbance at 280 nm)
<b>Storage</b>	For extended storage reconstitute product with 50% glycerol instead of water and then aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
<b>Note</b>	For research use only
<b>Application notes</b>	Not specified. Optimal titers for applications should be determined by the researcher.
<b>Expiration Date</b>	12 months from date of receipt.