

## Product Datasheet

# Fluorescein Conjugated Affinity Purified anti-Sheep IgG [H&L] [Rabbit] (orb1675069)

<b>Description</b>	Fluorescein Conjugated Affinity Purified anti-Sheep IgG [H&L] [Rabbit]
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Sheep
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FACS, FC
<b>Immunogen</b>	Sheep IgG whole molecule
<b>Preservatives</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Concentration</b>	1.0 mg/ml (by UV absorbance at 280 nm)
<b>Storage</b>	This product is stable for several weeks at 4° C as an undiluted liquid. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
<b>Note</b>	For research use only
<b>Application notes</b>	specified. Optimal titers for applications should be determined by the researcher.
<b>Expiration Date</b>	12 months from date of receipt.